EMS & TRAUMA SYSTEMS

Portland State Office Building 800 NE Oregon Street, Suite 305 Portland, OR 97232-2162





How to Make Hemorrhage Control Trainers for Stop The Bleed

- 1. Purchase a 6" x 36" foam roller (you can find them online approx. \$20.00). Do not buy the compressed Styrofoam rollers. These white rollers are softer and easier to create the manikin. They are also more lifelike.
- 2. Cut it in half—this will give you two training legs at 18" long each.
- 3. Cut or dig a jagged hole near one end of the leg that is approx. 2" in diameter and 3" deep.
- 4. Cover the leg in fabric (to simulate pants). Either sew fabric to make a tube or use an old pair of pants.
 - 5. Fold and hot glue the ends.
- 6. Cut a hole in the fabric over the wound to simulate a jagged gun or blast injury.
- 7. You can take a red Sharpie marker and color the wound (not shown here), but be advised, have students wear gloves, the red ink will rub off slightly and stain fingers.





