

Testing Guidance  
*Mycobacterium tuberculosis*

Oregon State Public Health Laboratory



TABLE OF CONTENTS

Green Headings: Instructions for Local Public Health Authorities

Purple Headings: Instructions for Private Care Laboratories

**SAMPLE SUBMISSION OVERVIEW** ..... 2

**AFB SMEAR AND CULTURE** ..... 2

    Local Public Health Authorities ..... 2

    Private Care Laboratories ..... 6

**AFB ISOLATE IDENTIFICATION** ..... 7

    Local Public Health Authorities ..... 7

    Private Care Laboratories ..... 7

**M. TUBERCULOSIS ISOLATE SUSCEPTIBILITY** ..... 9

    Local Public Health Authorities ..... 9

    Private Care Laboratories ..... 9

**GENEXPERT** ..... 11

    Local Public Health Authorities ..... 11

    Private Care Laboratories ..... 12

**REFERENCES** ..... 13

**APPENDIX A: AFB Smear/Culture – Additional Specimen Sources Accepted** ..... 14

**APPENDIX B: Reference Labs with Testing for Nontuberculous Mycobacteria**.... 15

## SAMPLE SUBMISSION OVERVIEW

The Oregon State Public Health Laboratory (OSPHL) will test appropriately submitted specimens from patients suspected of having tuberculosis. The OSPHL will test specimens submitted to help monitor therapy of previously diagnosed TB patients or to aid in contact investigations.

The OSPHL accepts primary specimens from Local Public Health Authorities (LPHAs) and their affiliates. Primary specimens from Private Care facilities will be accepted at the request of a State or Local Public Health Authorities. In addition to *M. tuberculosis* (MTB) isolates submitted to this laboratory, in compliance with OAR 333-018-0018, the OSPHL will accept unidentified AFB smear-positive isolates to rule out MTB if the referring laboratory cannot. Private sector submitters will be charged a fee for this service. Submitters should contact the OSPHL 503-693-4100 for additional information.

## AFB SMEAR AND CULTURE

### Local Public Health Authorities

Whenever possible, collect and ship specimens early in the week and transport for receipt at the OSPHL Monday through Friday morning. Specimens should be transported to OSPHL as soon as possible after collection. Avoid batch transport of a series of specimens. Holding specimens may compromise specimen integrity, delay testing due to overgrowth of non-tuberculous organisms and inhibit identification of MTB.

The most common specimens submitted for this order are expectorated sputum and urine. Additional specimen sources accepted are provided in Appendix A.

### Expectorated sputum

- All required collection and shipping materials are available from OSPHL.
- The best specimen is the first deep sputum obtained in the morning. If possible, observe the first sputum collected by the patient to ensure they understand how to collect and refrigerate the sputum. Resources for sputum collection are provided below. Patients should be instructed to:
  - Clear saliva from the mouth,
  - Rinse the mouth with water, and
  - Collect the sputum into the collection container.
- Collect a series of three to five single, early morning samples in the collection tubes provided. Use a separate tube for each sample; do not collect multiple samples into a single tube.
  - A volume of 5 to 10 mL is optimal. However, smaller volumes will be accepted. Results from samples of less than 1.0 mL will be reported with a qualifier.

- Sample volume greater than 11mL will be split, centrifuged, and then pooled for testing.
- Label the specimen with collection date and at least two unique patient identifiers (e.g., first and last name, date of birth, medical record number).
  - If multiple specimens are collected on the same day, the time of collection is also required on the specimen to properly link the specimen with the test order.
- Complete and submit a [General Microbiology Test Request Form](#) for each specimen.
- Securely cap tubes to avoid leakage.
- Decontaminate outside of specimen tube with a tuberculocidal disinfectant (if available) or 70% isopropyl alcohol.
- Refrigerate specimens at all phases of storage and transit, including at the patient's home, transit between the patient's home and the Local Public Health Authority, at the LPHA, and during transit to OSPHL.
- The following resources for sputum collection are provided in multiple languages:
  - Home Sputum Collection: A Step-by-Step Guide, offered by Public Health of Madison and Dane County, WI (videos):  
<https://www.youtube.com/playlist?list=PL09RoNDnObhOcYPpeYMTG3BpPHU5-Le6H>.
  - Instructions for Collecting Sputum Fact Sheet, offered by the Oregon Tuberculosis Program (handouts):  
<http://www.oregon.gov/oha/PH/DISEASES/CONDITIONS/COMMUNICABLEDISEASE/TUBERCULOSIS/Pages/factsheets.aspx>

## Urine

- Collect a series of three to five single, early morning, mid-stream urine samples. Do not collect multiple samples into a single tube.
- Aseptically transfer up to 35 mL of urine to provided tube.
- Secure lid against leakage.
- Label the specimen with collection date and at least two unique patient identifiers (e.g., first and last name, date of birth, medical record number).
  - If multiple specimens are collected on the same day, the time of collection is also required on the specimen to properly link the specimen with the test order.
- Complete and submit a [General Microbiology Test Request Form](#) for each specimen.
- Store refrigerated pending shipment.

## Specimen Shipping

- Ship samples according to current federal regulations. Guidance can be located on OSPHL's [Shipping and Transport webpage](#).
- Improperly packaged or leaking specimens present a hazard and may be rejected.
- OSPHL provides complete kits for packaging and shipping of tuberculosis test specimens.
- Label the specimen with collection date and at least two unique patient identifiers (e.g., first and last name, date of birth, medical record number).
  - Time of collection is required if multiple specimens were collected on the same day.
- Complete all sections of the [General Microbiology Test Request Form](#).
- If a different sterile container was used for collection, transfer specimen to a 50 mL Falcon tube.
- Secure the lid to prevent leakage during transport.
- Include absorbent material in the specimen transport bag.
- Place the tube sideways into the plastic specimen transport bag.
- Expel the air from the bag.
- Close bag by zipping the top seal.
- Put the fully completed General Microbiology Test Request Form into **outer** bag pocket.
- Refrigerate specimen pending shipment to OSPHL.
- Send specimens with fully frozen ice packs to maintain refrigerated temperatures during transport.
- Transport for receipt at the OSPHL as soon as possible after collection. Specimens should be received within 24 hours of collection but no more than 5 days after collection.

**GENERAL MICROBIOLOGY REQUEST**  
Oregon State Public Health Laboratory (OSPHL)  
7222 NE Evergreen Pkwy, Suite 110, Hillsboro, OR 97124  
Information: 503-693-4100

**PATIENT INFORMATION**  
Patient's name (last, first, middle initial): \_\_\_\_\_  
Date of birth (month/day/year): \_\_\_\_\_ Sex: ☐ Female ☐ Male Patient ID/Chart number: \_\_\_\_\_  
Race: ☐ American Indian or Alaska Native ☐ Asian ☐ Hispanic or Latino ☐ Black or African American ☐ White ☐ Not Reported or Unknown ☐ Native Hawaiian or Other Pacific Islander ☐ Unknown ☐ Declined  
Ethnicity: ☐ Not Reported or Unknown ☐ Declined  
Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ ZIP: \_\_\_\_\_  
County of Residence: \_\_\_\_\_ Outbreak number: \_\_\_\_\_  
Copy results to: ☐ County of Residence ☐ State Public Health ☐ Other Public Health

**PATIENT INSURANCE INFORMATION**  
Insurance carrier name: \_\_\_\_\_ Policy number: \_\_\_\_\_ Group ID: \_\_\_\_\_ Day/Week/CD (if used for test): \_\_\_\_\_

**SPECIMEN / ISOLATE INFORMATION**  
Date of collection: \_\_\_\_\_ Time of collection (for AM, PM): \_\_\_\_\_ Specimen submitted in (for modification): \_\_\_\_\_  
Original specimen source: ☐ Sputum ☐ Nasopharyngeal swab ☐ Sputum ☐ Tissue ☐ Urine ☐ Blood ☐ Other: \_\_\_\_\_ Is an isolate being submitted? ☐ Yes ☐ No

**Specimen Transport Bag Labels:**  
platinum pack SPECIMEN BAG  
REFRIGERATE  
BIOHAZARD PELIGROSO  
SPECIMEN ONLY

## Testing and Result Reporting

This order choice includes a primary smear, solid and liquid culture. GeneXpert, AFB PCR identification, and MTB isolate susceptibility testing will automatically be added when necessary.

- **Primary smear**
  - Determines presence or absence of AFB organisms in primary specimen.
  - Preliminary report issued within 2 working days of receipt.

- Initial positive results for new patients are called to the submitter on day of test result. Submitters will also receive a result released on the day of testing.
- Submitters should check Copia then contact the OSPHL if a report has not been received within 5 days of specimen receipt.
- **GeneXpert**
  - Determines presence or absence of MTB DNA in primary sputum samples, as well as mutations that confer resistance to Rifampin.
  - Cepheid GeneXpert® testing is FDA cleared for AFB smear positive and smear negative sputa from patients not on anti-tuberculosis therapy (TB-Rx) and who have a presentation consistent with active pulmonary tuberculosis.
  - Performed on the **first positive smear** from a previously undiagnosed patient with less than 72 hours of treatment for tuberculosis.
  - Approval may be obtained for testing smear negative specimens from patients that meet established criteria. Contact the OSPHL for additional information.
  - Preliminary report issued within 3 working days of receipt or approval to test, whichever is later.
  - Initial positive results for new patients are called and written reports released as results are available.
- **Solid media**
  - The OSPHL uses 7H10 solid media to grow and isolate AFB organisms.
  - A preliminary report is generated when acid fast bacilli are isolated. This report will include colony count and organism identification as MTB complex or *M. avium* complex based on detection of nucleic acids by the DNA PCR method or morphology consistent with a previous identification.
  - A final negative report of “No growth of AFB” is generated up to 9 weeks of incubation.
- **Liquid media**
  - The OSPHL uses MGIT liquid media to grow and isolate AFB organisms.
  - A preliminary report is generated when acid fast bacilli are detected by fluorescence microscopy. This report will include organism identification as MTB complex or *M. avium* complex based on detection of nucleic acids by the DNA PCR method or morphology consistent with a previous identification.
  - A final negative report of “No Growth of AFB” is generated after 6 weeks of incubation.

- **Susceptibility testing**
  - Antibiotic susceptibility testing is performed on **initial MTB culture positive isolates from new cases**.
  - Susceptibility testing is repeated at three-month intervals if a patient's samples continue to be reported as MTB culture-positive. Intervals are determined based on collection dates of specimens on which susceptibility testing was previously performed.
  - If susceptibility testing is already in progress for the same patient, the testing may not be performed on a specimen submitted.
  - "Resistant" susceptibility results are reported and called to the submitter as soon as possible. Results may be accessed via Copia. All written result reports are released within 24 hours of test completion.
  - All antibiotics are tested at the critical concentration for the test method in use (e.g., MGIT or agar proportion) as defined by the Clinical Laboratory Standards Institute (CLSI).
- **MGIT PZA**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of PZA at the critical concentration of 100 ug/mL
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTB complex (MTBC) identification.
- **MGIT IRE**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of Isoniazid (critical concentration 0.1 ug/mL), Ethambutol (critical concentration 5.0 ug/mL), and Rifampin (critical concentration 1.0 ug/mL)
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTBC identification.
- **Agar Proportion**
  - A solid agar-based susceptibility method that quantifies the proportion of resistant mutants within a population. When the percentage of resistant mutants is greater than 1% of the population, the isolate is categorized as "resistant."
  - Ofloxacin, a second line drug, is evaluated at the critical concentration of 2.0 ug/mL.
  - Results are typically available within four to five weeks of MTBC identification.

### Private Care Laboratories

Please contact the OSPHL (503) 693-4100 for approval and additional information regarding submission of AFB Smear/Culture specimens.

# AFB ISOLATE IDENTIFICATION

## Local Public Health Authorities

LPHAs should submit primary specimens for smear and culture as indicated above. LPHAs do not order this test.

## Private Care Laboratories

### Specimen Collection

- Primary specimens should be collected and processed according to the private laboratory's internal laboratory protocols and CLSI standards. AFB smear positive isolates may be submitted for the purpose of ruling in or out *M. tuberculosis* and *M. avium*.
- Submitters should retain subcultures of all submissions.
- Submit AFB culture isolates on Lowenstein-Jensen, 7H10, or other comparable solid AFB media. Failure to submit **pure cultures** may delay or adversely affect the identification process.
- AFB smear positive liquid media from automated TB systems, such as MGIT or BTA, will be accepted for testing.
  - To submit specimens from BTA systems, **do not** submit in BacT bottle.
  - Remove 5mL of material from the BacT bottle, spin for 15 minutes at 4150 rpm. Decant the supernatant, leaving approximately 0.5mL. Resuspend and transfer to a microcentrifuge tube, sealed with parafilm, for submission.
  - Specimens that cannot be spun may be submitted in conical tubes sealed with parafilm if needed.
  - If a liquid culture is submitted, evaluate DNA PCR results with caution. Results should be correlated with smear and culture morphology before releasing patient results to charting system.

### Specimen Shipping

- Ship samples according to current federal regulations. Guidance can be located on OSPHL's [Shipping and Transport webpage](#).
- Specimens should be packaged for secure transport
  - Parafilm the closure of the primary specimen container.
  - Place primary specimen container in secondary containment (biohazard bag, 15 mL conical tube, etc.)
  - Place appropriately contained specimen into a biohazard bag with absorbent material prior to shipping.



## Testing and Result Reporting

### ○ AFB PCR

- A fluorescent smear will be performed prior to testing to confirm the presence of AFB bacteria before testing is initiated. If the concentration of AFB bacteria is less than a standardized inoculum, the test will be reported as Unsatisfactory for testing.
- Specimens are generally tested on Wednesdays and Fridays.
- Indeterminate results are repeated once before reporting.
- Results will be called to submitters at the conclusion of testing.
- If *Mycobacterium tuberculosis* complex is not detected for AFB positive submissions, no further testing will be performed at the OSPHL.
  - If further identification will be performed by a reference laboratory, it is recommended that a new specimen be collected if possible. A list of laboratories that may provide further testing is provided at the end of this document.
  - If a new specimen cannot be collected, submitting laboratory should consider referral of the retained subculture in their laboratory.
  - If other options have been exhausted, the OSPHL may be able to forward the existing isolate or specimen to the submitting facility or provider at their expense.
  - If a facility that did not submit the specimen needs the specimen sent from the OSPHL to their facility, a patient authorization must be obtained following the OSPHL procedures. This should be the last option utilized whenever possible.

### ○ Susceptibility testing

- Antibiotic susceptibility testing is performed on **initial MTB culture positive isolates from new cases**.
- Susceptibility testing is repeated at three-month intervals if a patient's additional samples are MTB culture positive. Intervals are determined based on collection dates of specimens on which susceptibility testing was previously performed.
- If susceptibility testing is already in progress for the same patient, the testing may not be performed on an isolate submitted.
- "Resistant" susceptibility results are reported and called to the submitter as soon as possible. Results may be accessed via Copia. All written result reports are released within 24 hours of test completion.
- All antibiotics are tested at the critical concentration for the test method in use (e.g., MGIT or agar proportion) as defined by the Clinical Laboratory Standards Institute (CLSI).



- **MGIT PZA**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of PZA at the critical concentration of 100 ug/mL
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTB identification.
- **MGIT IRE**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of Isoniazid (critical concentration 0.1 ug/mL), Ethambutol (critical concentration 5.0 ug/mL), and Rifampin (critical concentration 1.0 ug/mL)
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTB identification.
- **Agar Proportion**
  - A solid agar-based susceptibility method that quantifies the proportion of resistant mutants within a population. When the percentage of resistant mutants is greater than 1% of the population, the isolate is categorized as “resistant.”
  - Ofloxacin, a second line drug, is evaluated at the critical concentration of 2.0 ug/mL.
  - Results are typically available four to five weeks of MTB identification.

## **M. TUBERCULOSIS ISOLATE SUSCEPTIBILITY**

### **Local Public Health Authorities**

LPHAs should submit primary specimens for smear and culture as indicated above.  
LPHAs do not order this test.

### **Private Care Laboratories**

#### **Specimen Collection**

- Primary specimens should be collected and processed according to the private laboratory’s internal laboratory protocols and CLSI standards.
- **Isolates must be identified as *M. tuberculosis* complex (MTBC) prior to submission to OSPHL for susceptibility testing.**
  - If an AFB isolate has not yet been identified as MTBC, refer to the “AFB Isolate Identification” guidelines above.
- Submitters should retain subcultures of all submissions.
- Submit AFB culture isolates on Lowenstein-Jensen, 7H10, or other comparable solid AFB media. Liquid culture that has been positively identified as containing *M. tuberculosis* DNA may also be submitted.

- Failure to submit **pure cultures** may delay or adversely affect the identification process.

### Specimen Shipping

- Ship samples according to current federal regulations. Guidance can be located on OSPHL's [Shipping and Transport webpage](#).
- Specimens should be packaged for secure transport
  - Parafilm the closure of the primary specimen container.
  - Place primary specimen container in secondary containment (biohazard bag, 15 mL conical tube, etc.)
  - Place appropriately contained specimen into a biohazard bag with absorbent material prior to shipping.

### Testing and Result Reporting

- **Susceptibility testing**
  - Susceptibility testing is repeated at three-month intervals if a patient's additional samples are MTB culture positive. Intervals are determined based on collection dates of specimens on which susceptibility testing was previously performed.
  - If susceptibility testing is already in progress for the same patient, the testing may not be performed on an isolate submitted.
  - "Resistant" susceptibility results are reported and called to the submitter as soon as possible. Results may be accessed via Copia. All written result reports are released within 24 hours of test completion.
  - All antibiotics are tested at the critical concentration for the test method in use (e.g., MGIT or agar proportion) as defined by the Clinical Laboratory Standards Institute (CLSI).
- **MGIT PZA**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of PZA at the critical concentration of 100 ug/mL
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTBC identification.
- **MGIT IRE**
  - A liquid based susceptibility method that evaluates growth of MTB in the presence of Isoniazid (critical concentration 0.1 ug/mL), Ethambutol (critical concentration 5.0 ug/mL), and Rifampin (critical concentration 1.0 ug/mL)
  - Resistant results will be confirmed by repeat testing.
  - Results are typically available within two to three weeks of MTBC identification.

- Agar Proportion
  - A solid agar-based susceptibility method that quantifies the proportion of resistant mutants within a population. When the percentage of resistant mutants is greater than 1% of the population, the isolate is categorized as “resistant.”
  - Ofloxacin, a second line drug, is evaluated at the critical concentration of 2.0 ug/mL.
  - Results are typically available four to five weeks of MTBC identification.

## GENEXPERT

### Local Public Health Authorities

#### Specimen Collection

- Pre-approval is required for testing smear negative samples. Please contact the OHA TB Program for approval at OSPHL for additional instructions 503-358-8516.
- GeneXpert is automatically performed on AFB Smear/Culture orders when the first smear is positive from a previously undiagnosed patient. This test may only be performed for patients who are:
  - Not on anti-tuberculosis therapy (TB-Rx) and who have a presentation consistent with active pulmonary tuberculosis, or
  - Have been on TB-Rx for less than 72 hours.
- All required collection and shipping materials are available from the OSPHL.
- The best specimen is the first deep sputum obtained in the morning. If possible, observe the first sputum collected by the patient to ensure they understand how to collect and refrigerate the sputum. Resources for sputum collection are provided below. Patients should be instructed to:
  - Clear saliva from the mouth,
  - Rinse the mouth with water, and
  - Collect the sputum into the collection container.
- Label the specimen with collection date and at least two unique patient identifiers (e.g., first and last name, date of birth, medical record number).
  - If multiple specimens are collected on the same day, the time of collection is also required on the specimen to properly link the specimen with the test order.
- Complete and submit a [General Microbiology Test Request Form](#).
- Securely cap tubes to avoid leakage.
- Decontaminate outside of specimen tube with appropriate a tuberculocidal disinfectant (if available) or 70% isopropyl alcohol.
- Refrigerate specimens at all phases of storage and transit, including at the patient’s home, transit between the patient’s home and the Local Public Health Authority, at the LPHA, and during transit to OSPHL.

- The following resources for sputum collection are provided in multiple languages:
  - Home Sputum Collection: A Step-by-Step Guide, offered by Public Health of Madison and Dane County, WI (videos):  
<https://www.youtube.com/playlist?list=PLO9RoNDnObhOcYPpeYMTG3BpPHU5-Le6H>.
  - Instructions for Collecting Sputum Fact Sheet, offered by the Oregon Tuberculosis Program (handouts):  
<http://www.oregon.gov/oha/PH/DISEASES/CONDITIONS/COMMUNICABLEDISEASE/TUBERCULOSIS/Pages/factsheets.aspx>

## Testing and Result Reporting

- GeneXpert
  - Determines presence or absence of MTB DNA in primary sputum samples, as well as mutations that confer resistance to Rifampin.
  - Cepheid GeneXpert® testing is FDA cleared for AFB smear positive and smear negative sputa from patients not on anti-tuberculosis therapy (TB-Rx) and who have a presentation consistent with active pulmonary tuberculosis.
  - Preliminary report issued within 3 working days of receipt or approval to test, whichever is later.
  - Initial positive results for new patients are called and written reports released as results are available.

## Private Care Laboratories

### Specimen Collection

- Private care facilities will be assessed a fee for this test. Please contact the OSPHL for billing information or if additional information or guidance is required, (503) 693-4100. OSPHL fees can be found at [www.healthoregon.org/phlbilling](http://www.healthoregon.org/phlbilling).
- Collect sputum according to internal laboratory protocol and CLSI guidelines.
- Process sputum using CDC NALC-NaOH decontamination methodology.

### Specimen Shipping

- A minimum of 1mL of the processed/concentrated pellet (stored and shipped at 2-8°C) should be submitted to the OSPHL for receipt within 2 days of the processing date. Ship samples according to current federal regulations. Guidance can be located on OSPHL's [Shipping and Transport webpage](#).

## Testing and Result Reporting

### ○ GeneXpert

- Determines presence or absence of MTB DNA in primary sputum samples, as well as mutations that confer resistance to Rifampin.
- Cepheid GeneXpert® testing is FDA cleared for AFB smear positive and smear negative sputa from patients not on anti-tuberculosis therapy (TB-Rx) and who have a presentation consistent with active pulmonary tuberculosis.
- Preliminary report issued within 3 working days of receipt.
- Initial positive results for new patients are called and written reports released as results are available.

## REFERENCES

1. U.S. Dept of Health, Education, and Welfare/Public Health Service/Centers for Disease Control. Public Health Mycobacteriology, A Guide for the Level III Laboratory. 1985.
2. U. S. Dept of Health, Education and Welfare/Public Health Service/Centers for Disease Control. Procedures for the Isolation and Identification of Mycobacteria. 1975. HEW Publication No. (CDC) 77-8230.
3. Centers for Disease Control and Prevention. 2009. United States Morbid. And Mortal. Weekly Rep. January 16, 2009 /58(01);7-10  
[http://www.cdc.gov/mmwr/preview/mmwrhtml/mm5801a3.htm?s\\_cid=mm5801a3\\_e](http://www.cdc.gov/mmwr/preview/mmwrhtml/mm5801a3.htm?s_cid=mm5801a3_e)
4. Centers for Disease Control and Prevention (CDC) GeneXpert MTB/RIF Assay Resource:  
<https://www.oregon.gov/oha/ph/DiseasesConditions/CommunicableDisease/Tuberculosis/Documents/tools/NAATGuide.pdf>
5. Official American Thoracic Society/Infectious Diseases Society of America/Centers for Disease Control and Prevention Clinical Practice Guidelines: Diagnosis of Tuberculosis in Adults and Children *Clin Infect Dis.* first published online December 8, 2016 doi:10.1093/cid/ciw694

## APPENDIX A: AFB Smear/Culture – Additional Specimen Sources Accepted

The most common specimens submitted for this order are expectorated sputum and urine. Additional specimen sources accepted are provided below.

- Gastric lavage
  - Specimens which cannot be processed within four hours of collection must be neutralized. For each 35 – 50 mL of specimen, use 1.5 mL of sterile 40% anhydrous disodium phosphate ( $\text{Na}_2\text{HPO}_4$ ) or two pH 7.4 pHydron buffer capsules or tablets.
  - Secure lid against leakage and label specimen tube with two patient identifiers, and collection date.
  - Complete and submit a General Microbiology Test Request Form for each specimen.
  - Store at refrigerated temperatures pending and during shipment.
- Body Fluids (cerebrospinal fluid, thoracentesis, pleural fluid, synovial fluid, etc.)
  - Aseptically transfer specimen to plastic specimen tube. There is no minimum volume.
  - Secure lid and label collection tube with two patient identifiers, and collection date.
  - Complete and submit a [General Microbiology Test Request Form](#) for each specimen.
  - Store at refrigerated temperatures pending and during shipment.
- Wounds and lesion
  - Aseptically obtain purulent material with a sterile swab or by aspiration or washing.
  - Aseptically transfer specimen to a sterile tube. Add a small amount of 7H9, 7H11 or sterile isotonic saline to prevent specimen dehydration.
  - Secure lid against leakage and label specimen tube with two patient identifiers, specimen source and collection date.
  - Complete and submit a [General Microbiology Test Request Form](#) for each specimen.
  - Store at refrigerated temperatures pending and during shipment.

## **APPENDIX B: Reference Labs with Testing for Nontuberculous Mycobacteria**

This list is provided for reference and information only to medical providers. The Oregon State Public Health Laboratory (OSPHL) and TB Program, Oregon Health Authority do not endorse the services of reference laboratories on this list.

If a new specimen can be collected or a laboratory can ship the retained subculture at their laboratory, please do so. If another specimen cannot be collected, the OSPHL may ship the specimen to the facility or provider that submitted the specimen. Payment of testing, specimen transport, and coordination of shipping specimens from OSPHL is the responsibility of the medical provider or facility, not local or state public health.

### **ARUP Laboratories**

Client Services: 800-522-2787

Lab Test Directory: <http://www.aruplab.com/testing>

### **Quest Diagnostics**

Customer Service: 866-697-8378

Test Center: <http://www.questdiagnostics.com/testcenter/TestCenterHome.action>

### **National Jewish Health Laboratories**

Main Customer Service: 200-550-6227, Option 6

Mycobacteriology (TB & NTM) Customer Service: 800-550-6227, Option 3

Test Catalog: <https://www.nationaljewish.org/professionals/clinical-services/diagnostics/adx/search-adx-tests>