

# ARIAS System Update

RELEASE NOTES

---



## Contents

**ARIAS Release Notes: 7.1.2020** .....2

    What's new:.....2

    Bug Fixes:.....2

## ARIAS Release Notes: 7.1.2020

### What's new:

- Added a new option for “Monitoring Status”. The new option “Marked for deduplication” is intended to be used when users run into duplicates after the import. Existing status options were found to not be a good fit for accurately describing these records.

The screenshot shows a user interface with a 'Details' section. Under 'Automated Monitoring', the value is 'Yes'. The 'Monitoring Status' dropdown menu is open, showing options: 'Monitoring', '--Select--', 'Completed Monitoring', 'Opted Out', 'Lost to Follow-up', 'Transferred to Orpheus', and 'Marked for deduplication'. The 'Marked for deduplication' option is highlighted in blue. Below the dropdown is a 'Last Assessment Date' field with a lock icon. At the bottom, the 'Address' section is partially visible.

- Date definitions in the Spanish assessment corrected from “(DD/MM/YYYY)” to “(DD/MM/AAAA)”.
- Assessment symptom “Fever (100 or above)” changed to “Fever” to align with new guidelines.
- Assessment temperature adjusted so any temperature value equal to or above “100.4” will automatically check the “Fever” symptom.

### Bug Fixes:

- Fixed: Issue was identified with some Sprint carrier phones splitting the automated text message into two separate messages. This caused the link to become invalid as the break occurred within the link. The link was moved earlier in the message, so it is maintained even if the message is split in these situations.