



CannaSafe

Independent Testing & Adulterants

D. Duncan, PhD¹ & P. Trujillo²

¹ Laboratory Director, CannaSafe

Tigard, OR

² Technical Director, ChemHistory

Milwaukie, OR

Outline

Rules & Testing

Terpene Analysis


Commercial Terpene Solutions

Adulterants in Inhalable Products

From Rules to Data: A laboratory prospective



CannaSafe



How to Quantify

Targeted Testing

- ▶ Low to moderate cost
- ▶ Chemically isolate target compounds from sample
- ▶ Easy to automate
- ▶ Multiple detectors available

Untargeted Testing

- ▶ Moderate to high cost
- ▶ Chemically isolate class of compounds from sample
- ▶ Difficult to automate
- ▶ Time of flight mass spectrometers





CannaSafe

Typical GC Chromatography




CannaSafe



Trends




CannaSafe



Typical Flavored Distillate Data




CannaSafe



Myrcene Terpene Solution




CannaSafe



Fruity Pebbles Terpene Solution



CannaSafe




Geraniol Terpene Solution




CHEMHISTORY
QUALITY CONTROL LABORATORY



CannaSafe



Geraniol Terpene Solution (zoom)



What have we seen?



CannaSafe

Sample	Squalene (%)	MCT (%)	Bisabolene (%)	Vitamin E (%)	Vitamin E Acetate (%)
1 Viscosity	<0.10	2.12	ND	ND	ND
2 Terp. Sol.	<0.10	7.31	ND	ND	0.47
3 OLCC RSO	<0.10	5.82	<0.50	<0.10	0.31
4 CBD Vape	<0.10	7.09	ND	<0.10	0.40
5 OLCC Vape	0.66	2.14	3.79	0.12	0.35
6 Viscosity	<0.10	ND	<0.50	ND	0.71
7 CBD Vape	<0.10	ND	0.88	<0.10	0.44

Summary

- ▶ Independent testing laboratories are focused on targeted testing
- ▶ Terpene analysis separates compounds by volatility
- ▶ Typical flavored distillate products have few detected compounds
- ▶ Low-volatility compounds are found in terpene solutions
- ▶ Non-cannabis terpenes were found in OLCC regulated product
- ▶ Many viscosity agents can be measured with current equipment



CannaSafe