

OREGON MILITARY DEPARTMENT	NUMBER: 99.200.13
ADJUTANT GENERAL PERSONNEL	EFFECTIVE DATE: Dec16, 2020
SUBJECT: Automatic External Defibrillator (AED)	

APPLICABILITY:

This policy and the procedures contained herein are applicable to state employees.

AUTHORITY/REFERENCE:

ORS 431A.455, ORS 30.802, ORS 339.345

ATTACHMENTS:

AED Quick Reference Guide.

PURPOSE:

The purpose of this policy is to provide information regarding the use of an AED.

GUIDANCE:

AED's will be reviewed during quarterly safety inspections ensuring that the AED is maintained in a state of readiness and maintenance is documented.

Maintaining Readiness:

Maintenance will be the responsibility of the ORARNG State Safety Office, with the exception of Portland Air National Guard Base (PANG), Kingsley Air National Guard Base, and Oregon Youth Challenge Program (OYCP). Maintenance at PANG and Kingsley Field is the responsibility of the Federal Medical Clinic on each base. Maintenance at OYCP will be the responsibility of a designated person at that facility. The AED will be checked for readiness after each use and at least once every 90 days if it has not been used in the preceding 90 days. Checks will include the following:

- Assure that the device is operational, fully charged battery and turns on (AED specific e.g., OK is visible in the readiness display).
- Check the expiration date on the electrode packet (AED specific). If the date has expired, replace.
- If used, take the AED out of use until it is checked, and pads are replaced.

Training:

- Employees are notified of the AED's location when hired and anytime the location changes.
- Employees are encouraged to find an online training in the use of AED's. There are several options for free online training such as

<https://www.firstaidforfree.com/online-automated-external-defibrillator-aed-course/>

INQUIRIES / QUESTIONS: Questions pertaining to this guidance may be directed to Adjutant General Personnel at (503) 584-3588.



Tracy Garcia
Director
Oregon Military Department
Adjutant General Personnel

Adult CPR Reference

1 Shake and shout



Call for help
local emergency number

Send someone
to get AED

2 Check for breathing




Open airway
Look, listen, feel



If not breathing,
give 2 slow breaths

3 Check signs of circulation
(pulse and/or normal breathing,
coughing or movement)

If no signs of circulation



Activate AED as soon as available
or
Perform CPR until AED is available




Compress and breathe at a ratio of 15 to 2

AED Reference

1 Turn on AED

2 Bare chest and apply electrode pads

3 Follow visual and voice prompts



Medtronic Physio-Control
11811 Willows Road Northeast
Post Office Box 97006
Redmond, WA 98073-9706
Tel: 1.800.442.1142
www.aedhelp.com

Europe
Tolochenaz, Switzerland
Tel: 41.21.803.8000
Fax: 41.21.803.5099

LIFEPAK® Family of AEDs



LIFEPAK CR Plus
defibrillator



LIFEPAK 500
automated external defibrillator



LIFEPAK 20
defibrillator/monitor



LIFEPAK 12
defibrillator/monitor

Model 500, CR Plus, and CR Plus are registered trademarks of Medtronic Physio-Control Corp.
©2006 Medtronic Physio-Control Corp. N/A 001128-000 / 1-06 20070706a